
Instructions for Use

Gebrauchsanweisung
Instrucciones de Uso

FlexLine

REF	010101
------------	--------


CONTENTS

ENGLISH 3
DEUTSCH 6
ESPAÑOL 9
SYMBOLS 11
INFO 12

ENGLISH

Read carefully before use

Keep for future reference

Components: 1. Driver nut, 2. Base, 3. Adhesive, 4. Inner guide

Intended use:

FlexLine is a needle guide intended to provide assistance in guiding and holding needles or electrodes during MRI, X-Ray or CT-guided needle-based interventions. The device enables physicians to set the needle or electrode in an angle and position according to the desired trajectory by the physician. FlexLine is able to keep the alignment of the needle or electrode during the intervention once the device is locked and adhered on the patient's skin.

Possible known complications:

Possible complications related directly with the medical device include minor irritations and allergic reactions due to the adhesive surface of the device.

Specifications:

- Inner guide's gauge sizes included with the device for 14G, 16G, 18G or 20G needles
- Possible angle settings in 2 degrees of freedom from -40° to 40°.
- 3M™ High Adhesion Double Coated Polyester Tape
- Material of the needle guide: High Impact Polystyrene (HIPS)
- Total height of the device: 40mm

Contraindications:

The use of FlexLine is contraindicated when/if:

- The skin's surface area where FlexLine is meant to be fixed is injured, has not been properly disinfected or is too humid.
- The packaging shows signs of wear or damage.
- Any component of FlexLine is loose or out of place according to the schematic figures.
- The needle or electrode's gauge size is not compatible with the inner guides provided with the device.
- The needle or electrode's weight exceeds the limit established in this document.
- FlexLine does not cover entirely its base's surface area on the patient's skin.
- The device has been used before.

Warnings:

06/2021 InLine-Med GmbH

- The medical device should only be used by qualified personnel in needle-based procedures, such as physicians or radiological technicians.
- Sterile device for one time use ONLY. Not to be reused or re-sterilized. The adhesiveness and sterility is not guaranteed after one use.
- The needle guide is intended to hold the alignment of needles or electrodes, but does not fix them. The needles and electrodes remain movable with regard to their depth.
- FlexLine is not intended to determine by itself the needle or electrode's entry point, this is determined by the physician.
- FlexLine does not provide any assessment on whether the selected entry point and alignment is correct or incorrect. This is purely the physician's decision.
- FlexLine does not provide any suggestions or recommendations to select the entry point and alignment. This is purely the physician's decision.
- The medical device is designed and manufactured to hold the weight of needles or electrodes. External additional forces or improper use of the device can lead to injuries or misalignment.
- Applying excessive force while locking the device may cause its breaking and can lead to injuries or misalignment.

Warning:

The company InLine-Med GmbH does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive can degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

The following steps must be performed prior to the application:

1. Check that the package has not been opened or damaged before use; check the sterility date.
2. Open the package.
3. Verify the size of the inner guide according to the needle's gauge size before use.

Directions for use:

1. Insert the inner guide (4) with the desired gauge size for the biopsy needle or electrode to be used.
2. Remove the adhesive paper (3) and place the device on the desired area for the intervention by placing the base (2) on the patient's skin.
3. Set the desired angle and lock the system through the driver nut (1) by rotating it in clock-wise direction.
4. Insert the biopsy needle or electrode to be held.

Disposal:

- Dispose the device as appropriate, immediately after use.
- Ensure that the device is disposed of safely and in accordance with all current national and international waste disposal directives.

Operating conditions:

Operate with a temperature between 5°C – 30 °C (41°F – 86 °F).
Humidity between 41% and 86%

Storage conditions:

Protect from moisture. Keep dry.
Store in a dark, dry place at a temperature between 5 – 30 °C (41°F – 86 °F).
Humidity between 41% and 86%

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